
## Quickstart Guide: Switch Click™ USB

ableNet®

www.ablenetinc.com

### Description

The Switch Click USB emulates a mouse or keyboard hit (default left mouse click) to a computer when activated. The device has a USB connection and switch input jack allowing you to attach your favorite switch to the Switch Click USB.

### Basic Operation:

1. Connect the Switch Click USB to any available USB port on your computer.
2. Press the switch top to send a left-click to your computer.

### Configuring Device Output:

The Switch Click USB can send one of several signals to a computer. To select the desired signal:

1. Ensure that the Switch Click USB is not plugged into a computer.
2. Remove the bottom panel by loosening the screws.
3. Configure the dip switches as shown at right.
4. Replace the bottom panel and plug the Switch Click USB into a computer.

Switch in ON position: .....	Mouse/Keyboard Function
None .....	Left Click (default)
1 .....	Right Click
2 .....	Space
3 .....	Enter (example at left)
4 .....	Tab
5 .....	Up Arrow
6 .....	Down Arrow
7 .....	Left Arrow
8 .....	Right Arrow

For best results use the Switch Click USB with scanning software. *Switch Click USB will not move a cursor or arrow around the screen.* Configure any scanning software to recognize the same signal as the Switch Click USB.

### Using an External Switch:

The switch jack on the side of the device allows you to connect any single switch (mono 1/8" plug). This jack performs the same function as the switch top. It does not add a second click behavior.

### Troubleshooting

If you have problems with the Switch Click, first ensure that any scanning software is correctly configured. Also consider opening your computer's control panel software, examining mouse functions.

### Warranty/Repair/Return Information

AbleNet offers a limited two-year warranty on the Switch Click product. For full warranty details, visit [www.ablenetinc.com](http://www.ablenetinc.com)

### For more information on AbleNet products and services:

AbleNet, Inc., | Minneapolis/St. Paul, MN | 800-322-0956 (US & Canada) | Internet: [www.ablenetinc.com](http://www.ablenetinc.com)  
© 2009 AbleNet, Inc., All Rights Reserved

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© 2009 AbleNet, Inc. Tous droits réservés
© 2009 AbleNet, Inc. Alle Rechte vorbehalten