PowerLink 4 control unit
User Manual

Read this entire Manual before operating this device.

English

Español

Français

Deutsch

Symbols

The exclamation point within an equilateral triangle is intended to alert the user to the presence of important operating and maintenance (servicing) instructions in the literature accompanying the appliance.

WARNING: indicates the possibility of injury associated with use or misuse of the device

CAUTION: indicates the possibility of a problem with the device associated with use or misuse.

Type BF Applied Part

WARNINGS

IMPORTANT SAFETY INSTRUCTIONS

WARNING: TO REDUCE THE RISK OF FIRE, ELECTRIC SHOCK OR PERSONAL INJURY:

- Turn off/unplug before cleaning or other maintenance
- The mains plug is used as the disconnect device, the disconnect device shall remain readily operable
- Make sure plug is readily accessible at all times
- Read these instructions
- Heed all warnings
- Follow all instructions
- Do not use this device near water.
- Install in accordance with manufacturer's instructions
- Do not operate the unit if the cord or unit has been damaged. Refer all servicing to qualified service personnel. Servicing is required when the apparatus has been damaged in any way, such as power-supply cord or plug is damaged, liquid has been spilled on or objects have fallen into the apparatus, the apparatus has been exposed to rain or moisture, the apparatus does not operate normally or it has been dropped.
- Do not operate the device in close proximity to flammable or combustible vapors, such as aerosol (spray) products or where oxygen is being administered.
- Do not block any ventilation openings.
- Do not defeat the safety purpose of the polarized or grounding-type plug. A polarized plug has two blades with one wider than the other. A grounding type plug has two blades and a third grounding prong. The wide blade or the third prong is provided for your safety. If the provided plug does not fit into your outlet, consult an electrician for replacement of the obsolete outlet.
- Protect the power cord from being walked on or pinched particularly at plugs, convenience receptacles, and the point where they exit from the apparatus.
- Only use attachments/accessories specified by the manufacturer.
- Unplug this apparatus during lightning storms or when unused for long periods of time.
- Apparatus shall not be exposed to dripping or splashing and no objects filled with liquids, such as vases, shall be placed on the apparatus.
- For indoor use only

CAUTIONS

- When in use, take caution to avoid pulling or creating tension on the cords.
- This device must be connected to a mains socket outlet with a protective earthing connection
- No user serviceable parts. Do not attempt to disassemble. This may cause the device to be damaged, malfunction or may cause personal injury
- This product can interfere with other electrical equipment if placed too close

SAVE THESE INSTRUCTIONS

Indications for Use

The PowerLink 4 is intended to provide switch access to/control of electrical appliances with basic on-off functionality. This product allows any one or two switches to turn an electrical appliance on or off for various durations.

Table of Contents


Product Map/Features	4
Using the PowerLink 4	
Basic Setup and Use	5
Linking with AbleNet Wireless Devices	5
Modes of Control	6
Direct	
Latched	
Timed	
2-Switch	
Count	
Troubleshooting	7
Maintenance	7-8
Care	
Cleaning	
Storage	
Technical Specifications	8
Regulatory	8
Warranty/Repair/Return	8

English

*Questions? Contact AbleNet, Inc. at 800-322-0956 or
e-mail us at customerservice@ablenetinc.com*

PowerLink Product Map & Features

Mode Button: Use this button to cycle between PowerLink's 6 Modes of Control. Mode LEDs provide visual feedback.

Status Buttons: Use these buttons for testing the PowerLink's connection to attached appliances. Status Buttons are also used when linking AbleNet Wireless devices. (See page 5.)

+/- Buttons: Use the +/- Buttons to program Timed Minutes or Timed Seconds Modes (see page 6).

Product Features

- Two grounded AC power inlets with flexible cords for customized positioning.
- AbleNet wireless compatibility with Jelly Beamer® and other AbleNet wireless devices.
- Two switch input jacks  (1/8 in/3.5 mm mono plug)
- Power ON/OFF button, Test buttons, 2-digit display
- 6 Modes of Control: Direct, Latch, Timed Minutes, Timed Seconds, Count, and Two-Switch. (See page 6 for details.)
- 120V models (US version) two outlets to allow choice-making activities with a variety of appliances, up to a total of 12 Amps.
- 230V models (Outside US version) two outlets to allow choice-making activities with a variety of appliances, up to a total of 10 Amps.
- No accessories provided.

Using the PowerLink 4

Basic Setup and Use

1. Plug the PowerLink into any AC wall outlet with the cord provided and turn the device ON.
2. Plug appliance(s) to either power cable on the side of the PowerLink and turn the appliance(s) ON.
Note: For pressure sensitive switches, consider tape or a Velcro® strap to hold the button down.
3. Test your AC power connection by pressing the appropriate Test button. The attached appliance should activate.
4. Attach any corded, single switch(es) to the corresponding switch jack(s) .
5. Press the Mode button until the "Direct" mode LED lights up.
6. Press and Hold switch(es) to activate their corresponding appliances. Release to turn OFF.
PowerLink's Status lights will activate and remain on whenever the PowerLink receives switch input.

Linking with AbleNet Wireless Devices

The PowerLink can be controlled by any AbleNet Wireless* transmitter such as the Jelly Beamer® wireless switch. Linked devices can access the PowerLink from up to 30 feet away without pointing or "line of sight".

Up to 8 wireless switches can be linked to one or both sides of the PowerLink, and will only activate sides that they are linked to.

To link wireless switches:

1. Turn both the PowerLink and wireless switch ON.
2. Press and Hold PowerLink's Status button on the side you wish to link for 6-10 seconds until the light turns solid Green.
3. Activate any AbleNet Wireless transmitter.
PowerLink's corresponding light will blink Red, then Green.
Repeat this step to link as many as 7 additional transmitters on that side.
4. Press and Release Status button to exit linking mode.
PowerLink's corresponding light will turn off.
5. Repeat steps 2-4 for opposite side.

NOTE: Once you exit linking mode, re-initiating the linking process will overwrite all previously linked devices.

To un-link wireless switches:

Press and Hold the desired Status for 10-15 seconds until the light blinks Red 3x, then Green, then off.

*AbleNet Wireless refers to devices using the AbleNet Wireless protocol and does not include Bluetooth™ devices such as the Impulse™ EMG switch, or infrared devices such as the AirLink® cordless switch. For a complete list of AbleNet Wireless switches, visit www.ablenetinc.com/wireless or contact AbleNet customer service.

Modes of Control

The PowerLink control unit offers six modes of appliance control: Direct, Count, Two-Switch, Timed Seconds, Timed Minutes, and Latch. Press the "Mode" button repeatedly to cycle through the modes of control.

Some modes use the Display Counter to convey information. These modes, and the display's function are identified with this symbol: .

Direct

The Direct mode keeps the appliance(s) on only while the switch is activated. Once released, the appliance(s) will turn off.

Count

Count mode behaves like Direct mode, showing the number of activations on the Display Counter.

 When used with two switches, the Display Counter will reflect whichever side was most recently activated. To view either side's number of activations, press the corresponding Status button.

Press and Release the Left and Right buttons simultaneously to reset the counter back to 00.

Two-Switch

Two-Switch mode requires two simultaneous switch activations to turn an appliance on. When activated, all attached appliances will activate in Direct mode.

Two switches must be attached or linked to use this mode.

Timed Seconds

The Timed Seconds mode allows the appliance to run for a preset amount of time, from 1 to 99 seconds, with a single activation.

Set the desired number of seconds with the Up and Down arrows. (Press and Hold to quickly move through the numbers.) The Display Counter will show the time set and, once a switch is activated, count down to one. Press the Up or Down arrow during activation to adjust the amount of time remaining.

 If a second switch is activated the Display will reflect that side's remaining time. Press either Status button to view the time remaining for that side.

Note: If you do not want to wait for the set time to expire, press the Mode button to select any other mode of control, attached appliances will shut off.

Timed Minutes

This mode behaves exactly like the Timed Seconds mode, but using minutes rather than seconds.

Latch

One activation of a switch turns the corresponding appliance on; the second activation turns it off. Use the Latch mode when you want an appliance (such as a lamp) to remain on indefinitely.

Troubleshooting

If the equipment fails to operate after you have read the instructions, follow the suggestions listed below. If the problem persists, please call AbleNet Customer Service at 1-800-322-0956 (US and Canada).

Do not attempt to dismantle this product. There are no serviceable parts.

Is PowerLink receiving power?

1. Does the wall outlet work with other devices?
2. Check to be sure the PowerLink is plugged into the wall outlet.
Check to be sure the power cord is plugged into the PowerLink.
3. Press the breaker button firmly to ensure that a circuit has not been broken. When broken, the breaker button will be extended, press this button back in until it clicks.
4. Turn on the PowerLink.

If the PowerLink lights turn on, move on to the next step.

Is the Appliance in good working condition?

1. Plug the appliance into a wall outlet.
2. Turn the appliance ON.
3. If the appliance has a pressure-sensitive switch (e.g., coffee grinder), tape the switch in the ON position.

If the appliance is working, re-attach it to the PowerLink and move on to the next step.

Does the PowerLink activate the Appliance with the Status buttons?

1. Put the PowerLink in "Direct" mode.
 2. Press the "Status" button on the side where the appliance is plugged into the PowerLink.
- If the appliance activates, move on to the next step.

Do switches activate the PowerLink and Appliance?

1. Attach any corded switch to the switch jack  on the same side as the Appliance.
2. Press the switch.

If the appliance activates, your PowerLink appears to be functioning normally. If you continue to have problems, contact AbleNet technical support at 800-322-0956.

Have any wireless switches/transmitters become un-linked from the PowerLink?

1. Check the wireless transmitter to ensure fresh batteries and power ON.
2. Make sure the wireless transmitter and the PowerLink are within 30 feet of each other.
3. Re-link wireless switches using the steps on page 5.

For Your Safety

1. Read the instructions for all appliances and devices before using them with the PowerLink.
2. Use caution when placing any appliance near a user (especially appliances that have moving parts or that generate heat, such as mixers and light bulbs).
3. Keep the floor clear of cords to avoid trips and falls, and to prevent the PowerLink or a connected appliance from being pulled to the floor. Position all cords near the wall. Always keep cords out of the reach of children.
4. Unplug the control unit when not in use and store it where it will not be pulled or dropped from a table accidentally.
5. Do not immerse the PowerLink control unit or any appliances in water. Do not place them in a location where they may get wet.
6. Do not open the PowerLink control unit to expose its interior.
7. To prevent injury, unplug the appliance(s) before opening the lid to a blender or food processor.
8. To prevent accidental activation of an appliance via a wireless transmitter or switch, unplug the PowerLink when not in use.

Maintenance

Care & Cleaning

Routinely check PowerLink before each use to confirm solid connection to AC power source. Do not use PowerLink if any cord appears to be damaged. If the device becomes wet or does not work properly, unplug the device immediately and do not attempt to use. Contact AbleNet, Inc. customer service at 800-322-0956.

- If cleaning is required, disconnect the PowerLink from the wall outlet; disconnect the appliances and switches from the PowerLink. The PowerLink should only be cleaned with a dry cleaning cloth. This device should never be immersed in liquid. This will damage the device electronics.
- Do not lift the PowerLink by the power cord.
- A PowerLink may be left connected to a frequently used appliance for the user's convenience; however, unplug the switch when not in use and note that Status buttons will continue to function and could cause accidental activation.
- Keep the device in a cool dry place. Keep the device away from heat, direct sunlight, and liquids. Do not store the device where it can fall or be pulled into water.

Technical Specifications

Electrical

Model: 100-10701
Electrical Class: Class I
Electrical Type: Class BF

Input Power/ Output Power:
Outside US- "230VAC 10A 50Hz"
Power Cord: US – Use only power cords rated 15A
Outside US – Use only power cords rated 10A

Disposal

Dispose of according to local regulations.

Operating Environment

Operating Temperature: 10 deg C to 40 deg C
Relative Humidity: 5% to 95% non-condensing

US – "120VAC 12A 60Hz"
Storage Temperature: 0 deg C to 40 deg C
Relative Humidity: 5% to 95% non-condensing

Transport and Storage

Switch Input

Only switches meeting the requirements of IEC 60601-1 shall be connected to the switch inputs.

Regulatory

EMC Compliance

This equipment has been tested and found to comply with Electromagnetic Compatibility (EMC) Requirement in accordance to EN60601-1-2. This device is compatible with its electromagnetic (EM) environment and does not emit levels of EM energy that cause electromagnetic interference (EMI) in other devices operating in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. Therefore, do not operate this device in an environment where other devices, which intentionally radiate electromagnetic energy in an unshielded manner, are being used. Portable or mobile radio frequency communication equipment can affect medical electrical equipment.

FCC Compliance

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Warranty/Repair/Return Information

AbleNet offers a limited two-year warranty on products. For full warranty details, visit www.ablenetinc.com.

For more information on AbleNet products and services:

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© 2009 AbleNet, Inc., All Rights Reserved
AbleNet, Inc., Minneapolis/St. Paul, MN
© 2009 AbleNet, Inc. Alle Rechte vorbehalten.