
Quickstart Guide: Switch Click™ USB

ablenet®

www.ablenetinc.com

Description

The Switch Click USB emulates a mouse or keyboard hit (default left mouse click) to a computer when activated. The device has a USB connection and switch input jack allowing you to attach your favorite switch to the Switch Click USB.

Basic Operation:

1. Connect the Switch Click USB to any available USB port on your computer.
2. Press the switch top to send a left-click to your computer.

Configuring Device Output:

The Switch Click USB can send one of several signals to a computer. To select the desired signal:

1. Ensure that the Switch Click USB is not plugged into a computer.
2. Remove the bottom panel by loosening the screws.
3. Configure the dip switches as shown at right.
4. Replace the bottom panel and plug the Switch Click USB into a computer.

| Switch in ON position: | Mouse/Keyboard Function |
|------------------------|-------------------------|
| None | Left Click (default) |
| 1 | Right Click |
| 2 | Space |
| 3 | Enter (example at left) |
| 4 | Tab |
| 5 | Up Arrow |
| 6 | Down Arrow |
| 7 | Left Arrow |
| 8 | Right Arrow |

For best results use the Switch Click USB with scanning software. *Switch Click USB will not move a cursor or arrow around the screen.* Configure any scanning software to recognize the same signal as the Switch Click USB.

Using an External Switch:

The switch jack on the side of the device allows you to connect any single switch (mono 1/8" plug). This jack performs the same function as the switch top. It does not add a second click behavior.

Troubleshooting

If you have problems with the Switch Click, first ensure that any scanning software is correctly configured. Also consider opening your computer's control panel software, examining mouse functions.

Warranty/Repair/Return Information

AbleNet offers a limited two-year warranty on the Switch Click product. For full warranty details, visit www.ablenetinc.com

For more information on AbleNet products and services:

AbleNet, Inc., | Minneapolis/St. Paul, MN | 800-322-0956 (US & Canada) | Internet: www.ablenetinc.com
© 2009 AbleNet, Inc., All Rights Reserved

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AbleNet, Inc., | Minneapolis/St. Paul, MN | 800-322-0956 (États-Unis et Canada) | Internet : www.ablenetinc.com
© 2009 AbleNet, Inc. Tous droits réservés
AbleNet, Inc., | Minneapolis/St. Paul, MN | 800-322-0956 (USA und Kanada) | Internet: www.ablenetinc.com
© 2009 AbleNet, Inc. Alle Rechte vorbehalten