

# Pilot One



Convert spoken commands into signals to control numerous devices such as TV, DVD, telephone, lights, door openers and more! The Pilot One features flexible voice programming, is portable, and offers easy programming after quick training from AbleNet. Pilot One can be used in any environment.

Note: The Pilot One is a sophisticated product involving substantial configuration. Contact AbleNet to locate a distributor near you for installation services.

## KEY FEATURES:

- 15 pre-set command menus including television, bed, TV, satellite, light and device
- Seamless operation of the Sero! Telephone and the IR Bed Controller
- No computer or laptop programming required
- Direct-to-unit Infrared training for easy programming
- New microphone settings for background noise adjustments just like the

Pair the Pilot One with these products for out-of-the-box functionality

### Sero! (Phone)

Item Number: 10011000  
Price: \$1195

### IR Bed Controller

Item Number: 82455  
Price: \$1135

## Pilot One

### Product Specifications:

#### Pilot One

---

##### Input

---

*Speech*

*Dual switch or single switch with scanning*

*Internal microphone*

*External keys for caregiver use*

---

##### Output

---

*Infrared transmitter for control of your household consumer electronics*

*Internal loudspeaker*

---

##### Auditory

---

*Voice feedback on menu commands*

**Item Number: 88460**

**Price: \$2100**

*Mounting systems available. Call for details.*

*For information on AbleNet's full line-up of Environmental Control Units visit [www.ablenetinc.com/ecu](http://www.ablenetinc.com/ecu)*



AbleNet, Inc.  
Roseville, MN USA

**Purchasing Info**  
Contact your AbleNet Sales  
Representative or order online at  
[www.ablenetinc.com](http://www.ablenetinc.com)

(800) 322-0956 [ablenetinc.com](http://ablenetinc.com)

# Pilot One



Convert spoken commands into signals to control numerous devices such as TV, DVD, telephone, lights, door openers and more! The Pilot One features flexible voice programming, is portable, and offers easy programming after quick training from AbleNet. Pilot One can be used in any environment.

Note: The Pilot One is a sophisticated product involving substantial configuration. Contact AbleNet to locate a distributor near you for installation services.

## KEY FEATURES:

- 15 pre-set command menus including television, bed, TV, satellite, light and device
- Seamless operation of the Sero! Telephone and the IR Bed Controller
- No computer or laptop programming required
- Direct-to-unit Infrared training for easy programming
- New microphone settings for background noise adjustments just like the

Pair the Pilot One with these products for out-of-the-box functionality

### Sero! (Phone)

Item Number: 10011000  
Price: \$1195

### IR Bed Controller

Item Number: 82455  
Price: \$1135