

Environmental Control Units

PowerLink® 4



(800) 322-0956
ablenetinc.com/powerlink



Get turned on!

The PowerLink® 4 control unit allows switch users to control up to two electrical appliances with single switches. Six unique modes of control allow you to define exactly how, and for how long, appliances will be turned on. PowerLink will accept any single switch, including AbleNet's Big or Jelly Beamer wireless switches.

PowerLink 4

Item Number: 100-10701
Price: \$229

KEY FEATURES:

- Easy set up with built-in display and a simple new design
- Powerful access with six modes of control
- Accurate programming and data collection like never before
- Flexible outlets allow custom positioning of up to two appliances, activating independently
- Linkable™ wireless technology for use with multiple Big or Jelly Beamer switches
- Two switch input jacks
- Great looks at a small size

MODES OF CONTROL:

- Direct: Behaves like a traditional switch.
- Count: Behaves like Direct, shows number of activations.
- Two-Switch: Behaves like Direct, both switches must be activated.
- Timed Seconds: Runs appliance from 1 to 99 seconds.
- Timed Minutes: Runs appliance from 1 to 99 minutes.
- Latch: First activation turns appliance ON; second activation turns it OFF.

Pair the PowerLink with these wireless switches for access from up to 30 feet away!



Jelly Beamer Transmitter

Item Number: 100-34200
Price: \$89



Big Beamer Transmitter

Item Number: 100-34300
Price: \$89