

Environmental Control Units

Primo!

by

possum

helping independence emerge



(800) 322-0956
ablenetinc.com/ecu



REMOTE CONTROL OF:

- TV, DVD, Satellite, and Cable boxes
- Lights and appliances*
- Telephone*
- Doors*
- Emergency contact - via phone, pager, alarm*
- Bed*
- Communication Aids*

Possum Primo! is a popular, easy-to-use, affordable dynamic device providing remote control of the immediate environment to give independence and reassurance for users and caregivers. Every Primo! system is tailored to meet the exact needs of each individual and can be adapted as users' needs change.

Primo! can be used as a simple system providing a limited number of functions or can be configured to control a comprehensive environmental control system.

Primo! is quick and simple to set up and enables easy user adjustments for volume brightness and scan speeds. Users may operate the Primo! by touching individual cells, any part of the screen, or external switch.



Stand sold separately

Primo!

Item Number: 700100900
Price: \$2700.00

DIMENSIONS & WEIGHT (APPROX)

Width 3.9 inches (100 mm)
Length 7.4 inches (190 mm)
Depth 1.7 inches (44 mm)
Weight 18.6 ounces (529 grams)

Primo! Stand

Item Number: 800100900
Price: \$79.00

**Accessory products required for control of lights, appliances, telephones, doors, etc. Contact us for further information.*

Primo! - Touchscreen Environmental Controller



Programmability

- Can be programmed to perform over 179 functions
- Easy to use, full on-board programming, no tools or hardware required
- Built-in infrared capture
- Macros (scenes) fully editable with up to 20 functions per macro
- Changeable telephone numbers and basic functions
- User-adjustable volume, brightness, and scan rate
- PC backup software with email option
- Custom template design, 18 "get started" device templates - Satellite, DVD, bed, etc.

Operation

- Touchscreen for direct access or use as a large switch for scanning
- Build up and send any satellite/teletext pages, digital TV
- Two auxiliary outputs for control of third party devices such as communication aids, alarms, computers, dispersed alarms, pager, etc.
- Two input jacks providing single or two-switch scanning access with auditory feedback
- Day or night mode

Display

- High Definition dynamic screen
- Repeater screen displays menu name and other useful information
- 14 clear symbols with text labels per screen

Phone Operation

- Internal phone book with 30 stored numbers with speed dial option
- Enter and dial any phone number
- Telephone symbols for mobile, home, and work numbers
- Macro build for home banking
- Last number redial & phone memory
- Compatible with infrared telephones including the Possum Sero!
- Programmable phone book via touch screen

Power and Warnings

- Powered via included wall adapter or rechargeable battery
- Up to 16 hours battery life (dependent upon usage) with fast recharge time
- Visual and audible low battery indication

Mounting Options

- Wall, floor, table, and wheelchair mountable
- Quick release mechanism for easy transfer



AbleNet, Inc.
Roseville, MN USA

Purchasing Info
Contact your AbleNet Sales
Representative or order online at
www.ablenetinc.com