



ablenetinc.com/brix



Communication Tools

Talking Brix™ Communicators



Easy as 1-2-3, still affordable at 4-5-6!

Simple communication has never been more affordable with AbleNet's Talking Brix! Use one Brix for personal reminders, or attach as many as you like, creating simple, scalable communication grids. Talking Brix are thin, light communicators with built-in magnets, perfect for carrying in a pocket, or placing around the room! For table-top users, Talking Brix use an ingenious tab and slot connector to link to other Brix. Create multi-message communicators in any arrangement you like. At a cost-conscious price, one, two, three or more messages are within your budget!

KEY FEATURES:

- Connectable devices for customizable layouts
- Easy single-message recording on each Brix
- 10 seconds of recording time
- Power on/off switch
- Rechargeable battery (built-in)
- 1.8" Activation area
- Free Snap Switch Cap included

Talking Brix

Dimensions

2.56" x 2.56" x 0.6" (L x W x H)

Activation Area

1.8"

Item Number: 100-03400



Purchasing Info

Contact your AbleNet Sales Representative or order online at www.ablenetinc.com/brix

