

Communication Devices



(800) 322-0956

www.ablenetinc.com/communication



AbleNet's AAC products are used by millions of people worldwide with a variety of medical and developmental conditions. From children with developmental language disorders, to individuals with severe physical disabilities that limit speech, AbleNet communication devices are available in a range of sizes and features.

When choosing a device, consider how many messages will be required, what access methods will be used, and whether the device is for a single person, or a group environment. For example, a strategically placed BIGmack® may be used by multiple users in a group home for reliably identifying shared needs. Meanwhile, a SuperTalker™ will grow with a child as his or her communication needs change over time.

For additional information or consultation, please contact AbleNet at (800) 322-0956.



Video Brix™



Talking Brix™



BIGmack® & LITTLEmack®



TalkTrac



BIG Step-by-Step™ & LITTLE Step-by-Step™



Step-by-Step™ with Levels & Step-by-Step™ Choice with Levels



Lingo™ All-in-One-Kit



iTalk2™ & iTalk2™ with Levels



QuickTalker™ Communicators (7, 12, 23)



SuperTalker™ Progressive Communicator



SoundingBoard™ v.3.1.1



Jive!™ Handheld Communicator

MESSAGES : LEVELS	1 : 1	1 : 0	1 : 0	4:2	Unlimited Sequential : 0	Unlimited Sequential : 3	7 single, 1 sequential : 0	2 : 0 (2 : 3)	5	1, 2, 4, or 8 : 8	3 x 3 per board	Up to 100 equal cells per page / up to 10,000 pages
RECORDING TIME	90 sec.	10 sec.	2 minutes	80 sec.	2 minutes	4 minutes	7 @ 25 sec., 1 @ 60 sec.	2 minutes	6 min. / 12 min. / 20 min.	16 minutes	Unlimited	Up to 14 hours
SWITCH INPUT	No	No	1	N/A	1	1	No	2	No	8	Blue2 or Bluetooth keyboard	1 USB
TOY / APPLIANCE OUTPUT / ECU	No	No	1	N/A	1	1	No	2	No	2	No	IR Capability
ACTIVATION AREA	165 in. x .54 in. – .42 cm x 1.37	1.8 inches	5 inches / 2.5 inches	.5 inches / .7 inches	5 inches / 2.5 inches	2.5 inches	0.6 inch	2.5 inches each	Static Locations: 2.167x1.20 in. QT7: 3.38x2.65 / QT12: 2.17x1.68 / QT23: 1.20x1.20 in.	1.75 inches up to 7.75 inches	Varies with board and device	3.5 inches x 2 inches
ACTIVATION PRESSURE	16 oz. – 453.59g	2 oz.	3-5 oz. / 2-3 oz.	5 oz.	3-5 oz. / 2-3 oz.	2-3 oz.	5-10 oz.	2-3 oz.		8 oz.	Finger tap	1 oz.
WEIGHT WITH BATTERY	1.59 oz. – 45 g.	1.5 Oz.	15.5 oz. / 8.5 oz.	.120 oz.	15.5 oz. / 8.5 oz.	8.5 oz.	4.5 oz.	12.5 oz.		32 oz. (2 lbs.)	Varies with device	12 oz.
BATTERY TYPE	Rechargeable	Rechargeable	9-volt	Rechargeable NiMH	9-volt	9-volt	Rechargeable Lithium Ion	9-volt	(4) AA Alkaline Batteries	4 AA	Rechargeable battery (built-in)	Lithium-ion (8 hrs. of use)
ACCESSORIES AVAILABLE	No	No	Travel Strap, Mounting Plate, Additional switch top colors (business casual or semi-formal)	N/A	Travel Strap, Mounting Plate, Additional switch top colors (business casual or semi-formal)	Travel Strap, Mounting Plate, Additional switch top colors (business casual or semi-formal)	Includes: 2-4-8 keyguards, Belt Clip, Travel Strap, Wrist Mount	Travel Strap, Mounting Plate	No	Travel Strap, Mounting Plate, Adjustable Stand	No	Jacket including 2 white stylus / Carrying Case
SPECIAL FEATURES	provide instant access to communication anywhere	Symbol Cover	Snap Switch Caps, 4 colors of Twist-off switch tops included	Charger included, Can use standard wristbands	Snap Switch Caps, 4 colors of Twist-off switch tops included	Snap Switch Caps, 4 colors of Twist-off switch tops included	Ultra-portable, Headphone Jack	Built-in Symbol Storage	Built-in symbols storage Flexible recording time by location	Built-in Overlay Storage and Step-by-Step functionality	Use the built-in AbleNet symbol library for instant access to hundreds of symbols, Support for auditory scanning for users with low vision	Scrolling Cells / Character, Word, and Sentence prediction / IDV4Life Vocab library, Includes Picture Symbol Set
MOUNTING CONNECTION (PRICE)	Internal Magnet	Internal Magnet	100-LTMP (\$22.00) 100-STMP (\$22.00)	N/A	100-LTMP (\$22.00) 100-STMP (\$22.00)	100-STMP (\$22.00)	N/A	100-STMP (\$22.00)	No	100-40300 (\$185.00)	iDevice Mounting Plates	N/A
HCPC CODE	E2500	E2500	E2500	E2500	E2500	E2500	E2500	E2500	E2500 (7), E2502 (12), E2502 (23)	E2502	-	E2510
ITEM NUMBER	10003504	10003400	10002100 (BIG) 10002400 (LITTLE)	10002000	10002200 (BIG) 10002500 (LITTLE)	10002600 10003600	10003200	10002700 (iTALK2) 10003300 (w/Levels)	10003501 (7), 10003502 (12), 10003503 (23)	10002800	available on the App Store purchase in iTunes	10011200
PRICE	\$59.00	\$129.00	\$119.00	\$129.00	\$165.00	\$189.00	\$199.00	\$149.00/\$199.00	\$159/\$179.00/\$199.00	\$359.00	\$49.99	\$5900.00