



Environmental Control Unit (ECU) Selection Guide

Go to www.ablenetinc.com for complete line of ECU Devices or call 800-322-0956

AbleNet offers a broad array of environmental control options to meet the wide-ranging needs of persons with disabilities. To choose the right ECU for you, consider the key features identified below. If you have questions about the features or use of any product, contact AbleNet. We'll ensure that you're comfortable with your decision, and guide you to a distributor in your area.



PowerLink® 4



Relax Series



Primo!™



Jive!™



Pilot™

Easy Tech

Touch Screens

Voice-Activated

	Easy Tech		Touch Screens		Voice-Activated
INPUT	Switch (Wired or Wireless)	Switch (Multiple or Scanning)	Direct, Scanning	Direct, USB Switch, Scanning	Voice, Switches (with adapter)
PRIMARY OUTPUT	AC Control	Infrared or Infrared + RF (X10)	Infrared	Infrared	Infrared or Infrared + RF (X10)
OTHER OUTPUT	—	X10 with Converter	INSTEON with Converter	INSTEON with Converter	INSTEON or X10 with Converter (depending on model)
COMPATIBLE DEVICES	AC Controlled Appliances	TV, DVD, Stereo, AC Appliances Etc.	Sero! Telephone, IR Bed Controller, Any IR Device, AC Appliances	Sero! Telephone, IR Bed Controller, Any IR Device	Sero! Telephone, IR Bed Controller, Any IR Device, AC Appliances
NUMBER OF DEVICES IT CAN CONTROL	2	1, 4 + X10	Many	Many	15 to Many
MOUNTING	NA	66160	Included 80010900 (Optional Stand)	NA	Included
FEEDBACK	Visual	Auditory, Visual	Auditory, Visual	Auditory, Visual	Auditory, Visual
ITEM #	10010701	82050, 82060, 82070, 82150, 82000, 82100	700100900	10011200	88405 88460
PRICE	\$229	\$369 - \$893	\$2700	\$3540	PRO - \$5950 ONE - \$2100