

Infrared Head-Operated Mouse Comparison



Product	TrackerPro	HeadMouse Extreme	SmartNav 3
Manufacturer	Madentec	Origin	NaturalPoint
<i>Plug & Play</i>	√	√	<i>No (Software Required)</i>
<i>Windows</i>	√	√	√
<i>Mac OS X</i>	√	√	<i>(partial)</i>
<i>AAC Devices</i>	√	√	–
<i>Switch Inputs</i>	√ (2)	√ (2)	√ (2)
<i>Splitter Cable</i>	√	–	√
<i>High-Impact Design</i>	√	√	–
<i>Adjustable</i>	√	–	√
<i>Size</i>	<i>1.75" DIA</i>	<i>3.7" x 2.2" x 0.5"</i>	
<i>Robust in Sunlight</i>	√	–	–
<i>Vertical Mount</i>	√	√	–
<i>Horizontal Mount</i>	√	–	√
<i>OEM'd into other products</i>	√	–	–
<i>Pin-point Accuracy</i>	√	√	–
<i>Price</i>	\$995	\$995	\$499
<i>Warranty</i>	<i>2 Years</i>	<i>2 Years</i>	<i>1 Years</i>

