



TrackerPro Overview

Head-operated Mouse



What is TrackerPro?

- A mouse operated by head movement
- A camera mounts on your computer & tracks a small reflective dot on your forehead or glasses



How is it installed?

- Simply plug it into your computer's USB port - just like a mouse
- No software installation necessary
- Affix a dot to your forehead or glasses
- Plug and play!



What are its key features?

- Hands-free cursor control
- Pin-point accuracy
- Dual switch input
- Ultra-rugged (solid aluminum)
- Very flexible mounting
- Adjustable (up and down)
- Universal (any device with USB)
- Smooth tracking
- Small & attractive
- No software required
- Sensitivity control
- Reliable



Mounting

- The most flexible mounting of any head tracking device
- Vertical, horizontal
- Desktop, laptop, tabletop



Who Needs It?



- Anyone who can't use their hands to move a mouse (but have good head control)
- Typical types of disabilities:
 - Spinal cord injury
 - Multiple Sclerosis
 - Mild CP
 - ALS



Clicking Solutions

- Any switch (or two)
- WISP (sip/puff switch)
- Magic Cursor (dwell software)



Care and Maintenance

- Protect the camera lens
- Clean as you would any lens
- Keep the reflective dot fresh



Sensitivity Adjustment

- “Turbo” switch
- Up is faster
- Down is slower
- Finer mouse speed adjustments can also be made in the mouse control panel
- Typical settings are Up on a Mac computer, down for Windows



What's in the box?



More Information

Purchase on-line at

www.ablenetinc.com

Or

www.madentec.com

