



(800) 322-0956
 ablenetinc.com/ecu



Environmental Control Units

Jive!

by

possum

helping independence emerge



Jive! is a totally new and proprietary handheld communication aid. Jive! offers a wide range of symbol and text communication, uses the incredibly clear Acapela voice synthesizer, and introduces full infra-red environmental control for the first time in an AAC (Augmentative & Alternative Communication) device of this size. Set up the Jive! either onboard or via the software included with purchase.

KEY FEATURES:

- Handheld, lightweight and portable
- Long battery life
- Touch accessible, highly visible screen - even in bright sunlight
- Comprehensive message building capabilities
- PCS Symbols
- Character, Word, and Sentence prediction
- 12 clear, synthesized, high quality Acapela voices
- Digitized (recorded) speech and sounds
- Programmable - on Jive! or via included PC software
- IDV4Life vocabularies

COMMUNICATION VIA:

- Single cell, automatically spoken messages
- Message building via symbols or text, across multiple cells and pages
- Learns word and sentence patterns for quick message generation
- Abbreviation & alternative expansion

INFRA-RED FROM CELLS ON PAGES

- Synthesized code library such as the Primo! - with meaningful words associated to codes.
- Learning Infra-Red
- Assignable favorite devices
- Scenes including page changes and sleep delays
- Multiple parameters available for Infra-Red

Jive! Package:

- Possum Jive! Communication Aid
- Power adapter
- 2 Styluses
- USB cable
- PC setup Software and User Manual
- Quick Start Guide

Item Number: 10011200
 Price: \$5900.00

DIMENSIONS & WEIGHT (APPROX)

Width 3.3 inches (85 mm)
 Length 5.7 inches (145 mm)
 Depth 1.5 inches (40 mm)
 Weight 12 ounces (350 grams)

Jive! – Handheld communicator with Environmental Control




Information Bar
Messages are built here. Selecting the info bar will cause the Jive! to speak.

Command Bar
Universal commands are accessed from the command bar. The command bar is fully editable.

Layout
A Layout is a collection of interlinked pages containing a user's vocabulary. Cells can navigate to pages, build a message, transmit infra-red and a lot more besides.

Scrolling Cells

Introducing the new “scrolling cells” feature, unique to the Possum Jive! Scrolling cells add a new dimension to AAC. Multiple cells are stacked on top of each other minimizing page navigation. The selectable scrolling cells function allows the user to easily access a variety of phrases for speedy message building, related messages on a topic, variations in verb tense and much more!



Message Building

- Character Prediction - helps you spell
- Word completion - helps you finish words you've started typing
- Context Based Word Prediction - suggests the next word you're going to use - based on context and usage
- Sentence Prediction - suggests sentences you've already used, sorted by recency.
- User added words and sentences are fully editable
- Abbreviation Expansion - use a few characters to generate a message
- Alternative Replacement - decline irregular verbs transparently with symbols
- Built-In calculator functionality

IDV4Life

- Coherent set of vocabularies to help develop communication and access
- 5 vocabularies - Logically progressing from simple system usage all the way to complex message building
- Suitable for children and adults
- Incorporates LifeSkills concept and develops it through to full home control
- Coherent color key
- Consistent look and feel throughout the layouts



Acapela Text-To-Speech from Acapela Group.
The Picture Communication Symbols © 1981-2008 by Mayer-Johnson LLC. All Rights Reserved Worldwide.
Used with permission.
Powered by Keypoint Technologies (UK) Limited.



AbleNet, Inc.
Roseville, MN USA

Purchasing Info

Contact your AbleNet Sales Representative or order online at www.ablenetinc.com