



# Blue2

*Bluetooth switch*

---



The *NEW* Blue2 Bluetooth switch provides single or dual switch access to iDevices running iOS 7 or newer, Apple desktop or laptop computers running OS X Mavericks or newer, and all switch accessible apps or software on iOS, OS X, Windows, and Android! Blue2 connects to your favorite device via a Bluetooth connection that is easy and quick to set up. New features include an improved switch activation surface that takes 40-percent less pressure to activate, multiple modes of operation for compatibility to all apps, a long-lasting rechargeable battery, and a sleek, modern design

**Item:** 10000017  
**MSRP:** \$179



**AbleNet, Inc.**  
www.ablenetinc.com  
(800) 322-0956

# Blue2

## Bluetooth switch

---



### PRODUCT FEATURES

- › Bluetooth wireless connectivity
- › Two built in switch activation areas
- › Two external switch jacks for use with almost any wired switch
- › Three pre-programmed modes of operation
- › Three programmable modes of operation - record one to four keystrokes per switch activation
- › Rechargeable battery

### Compatible with:

- › iOS Switch Control and most switch accessible apps
- › Mac OS X Switch Control and most switch accessible software
- › Windows switch accessible software
- › Android switch access and switch accessible apps
- › Other software that uses keystrokes for switch clicks